

LEPU MEDICAL




Disposable Endocutters & Disposable Endocutter Reloads

JQIIA, JQIIB, JQIIL

Outstanding Performance
Effective Operation

CE
0197


EN ISO
13485: 2016



1. Connector 2. Joint Arrow 3. Lock Button 4. 360°Rotation Knob; 5. Safety Button ;
6. Reset Button; 7. Handle; 8. Trigger Handle; 9. Rotation Knob

*Disposable endocutter has three length types (shaft+handle), A type (standard endocutter)=370mm, B type (short endocutter)=320mm, L type (long endocutter)=470mm.

JQIIA JQIIB JQIIL



CE 0197 EN ISO 13485: 2016

Catalogue No.	Cartridge Color	Total Length	Suturing Length	Suturing Thickness			Accessory type
				first row	second row	Third row	
JJIID301	Brown	202	30	2.0	2.5	3.0	Rotating type
JJIID302	purple	202	30	3.0	3.5	4.0	Rotating type
JJIID303	black	202	30	4.0	4.5	5.0	Rotating type
JJIID451	Brown	217	45	2.0	2.5	3.0	Rotating type
JJIID452	purple	217	45	3.0	3.5	4.0	Rotating type
JJIID453	black	217	45	4.0	4.5	5.0	Rotating type
JJIID601	Brown	232	60	2.0	2.5	3.0	Rotating type
JJIID602	purple	232	60	3.0	3.5	4.0	Rotating type
JJIID603	black	232	60	4.0	4.5	5.0	Rotating type

Jiangsu Brightness Medical Devices Co., Ltd.

Factory: No.66 Hehuan Road, Zhonglou Economic Development Zone,
Changzhou, Jiangsu Province, P.R.China

Shenzhen Office: No.9 Building, BaiWangXin Hi-Tech Park, SongBai Road,
XiLi, NanShan District, Shenzhen, China 518055