



Hygeia Medical

Ireland Germany Malaysia Indonesia
www.hygeia.icu



HYGEIA MEDICAL LIMITED

Kyletologue House, 9 Cherrywood
Dundalk 21222, Co. Louth, Ireland

HYGEIA MEDIZINTECHNIK

PF 2513, Hindenburgstrasse
46395 Bocholt, Germany





HYGEIA MEDICAL
INSTRUMENTS

Innovative Solutions In Medical Devices



Section 1

Laryngoscope Blades

Laryngoscope Handles

Laryngoscope Packing



Code :STD 101

LED BULB



Specification

The Flex tip is a Macintosh laryngoscope blade with a tip which is adjustable through 70° A lever controls the tip angle during incubation to lift the epiglottis, giving a clear view of the vocal cords. Contact with upper row or teeth and unnecessary pressure on the tongue and soft tissue are avoided.

SMEGA Hi-LIGHT fiber optic bundle: 4mm dia

SMEGA Xenon Halogen Technology: For up to 40% more light than conventional halogen bulb

SMEGA LED Bulb: High light LED bulb with white light life 5 years don't need for replacement

- SMEGA Hi-light Fiber optic bundle with 4.3mm dia. And up to 6500 individual micro fiber for improved light transmission and longer life.
- SMEGA LED Technology: For up to 40% more light than conventional halogen bulb
- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- These blades need no maintenance ever
- Both ends of the fiber bundle are protected for life
- 2.5V battery power for even more light
- Compatible with all handles to the ISO 7376 (green Standard)
- Autoclavable.
- Fits all standard F.O. handles (green standard)

Size Measurement

Code	Blade	B/No	Size
STD:101-2	Child	2	102 MM
STD:101-3	Medium Adult	3	130 MM
STD:101-4	Large Adult	4	155 MM

Require: 2 x C size batteries Supply set with hard case packing.



SMEGA Flex Tip + Fiber Optic (F.O) Laryngoscope Blade



LED BULB

Available Blads With Conventional & Mega C.O. Style

Code :STD 102

Specification

The Flex tip is a Macintosh laryngoscope blade with a tip which is adjustable through 70° A lever controls the tip angle during incubation to lift the epiglottis, giving a clear view of the vocal cords. Contact with upper row of teeth and unnecessary pressure on the tongue and soft tissue are avoided.

SMEGA Hi-LIGHT fiber optic bundle: 4mm dia

SMEGA Xenon Halogen Technology: For up to 40% more light than conventional halogen bulb

SMEGA LED Bulb: High light LED bulb with white light life 5 years don't need to be replaced

- :- Swivelling blade tip facilitates difficult incubation.
- :- Reduced likelihood of damage to teeth and soft tissue.
- :- Integrated Fiber Optics for the best possible illumination.
- :- High Voltage stainless steel construction.
- :- Autoclavable.
- :- Fits all standard F.O. handles (green standard)

Require: 2 x C Size batteries

Supply set with hard case or Packing.

Size Measurement

Code	Blade	B/No	Size
STD:102-2	Child	2	105 MM
STD:102-3	Medium Adult	3	130 MM
STD:102-4	Large Adult	4	155 MM

Require: 2 x C size batteries Supply set with hard case packing.





LED BULB

Code :STD 103



Size Measurement

Code	Blade	B/No	Size
STD-103-00	Premature	00	68 MM
STD:103-0	Neonate	0	80 MM
STD:103-1	Infant	1	92 MM
STD:103-2	Child	2	115 MM
STD:103-3	Medium Adult	3	135 MM
STD:103-4	Large Adult	4	155 MM
STD:103-5	x-Large	5	175 MM

Specification

- :- SMEGA Hi-light Fiber optic bundle with 4.3mm dia. And up to 6500 individual micro fiber for improved light transmission and longer life.
- :- SMEGA LED Technology: For up to 40% more light than conventional halogen bulb
- :- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- :- These blades need no maintenance · ever
- :- Both ends of the fiber bundle are protected · for life
- :- 2.5V battery power for even more light
- :- Compatible with all handles to the ISO 7376 (green Standard)
- :- Autoclavable.
- :- Fits all standard F.O. handles (green standard)

Require: 2 x C size batteries Supply set with hard case packing.






Size Measurement

Code	Blade	B/No	Size
STD:104-0	Neonate	0	80 MM
STD:104-1	Infant	1	92 MM
STD:104-2	Child	2	115 MM
STD:104-3	Medium Adult	3	135 MM
STD:104-4	Large Adult	4	155 MM
STD:104-5	x-Large	5	175 MM

Specification

- SMEGA Hi-light Fiber optic bundle with 4.3mm dia. And up to 6500 individual micro fiber for improved light transmission and longer life.
- SMEGA LED Technology: For up to 40% more light than conventional halogen bulb
- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- These blades need no maintenance · ever
- Both ends of the fiber bundle are protected · for life
- 2.5V battery power for even more light
- Compatible with all handles to the ISO 7376 (green Standard)
- Autoclavable.
- Fits all standard F.O. handles (green standard)





LED BULB

Code :STD 105



Specification

The Miller FixTube Fiber Optic (F.O) laryngoscope blades with improved geometry simplify intubation. The shortened shoulder or the blade dramatically improves the anesthetist's view of the epiglottis and vocal cords. The distal lip or the blade features a special atraumatic shape

SMEGA Hi-LIGHT fiber optic bundle: 4mm dia
SMEGA LED Technology For up to 40% more light than conventional halogen bulb

- Smooth shape avoid, sharp edges and cavities for easy cleaning and sterilizing.
- These blades need no maintenance ever
- Both ends of the fiber bundle are protected for life
- 2.5v battery power for even more light
- Compatible with all handles the ISO 7376 (green Standard)
- Autoclavable.
- Fits all standard F.O handles (green standard)

Size Measurement

Code	Blade	B/No	Size
STD:105-00	Premature	00	65 MM
STD:105-0	Neonate	0	78 MM
STD:105-1	Infant	1	102 MM
STD:105-2	Child	2	154 MM
STD:105-3	Medium Adult	3	195 MM
STD:105-4	Large Adult	4	205 MM

Require: 2 x C size batteries Supply set with hard case packing.





LED BULB

Code : STD:106



Size Measurement

Code	Standard Set Details
STD-106-11	3 Blades Set 00 ,0 ,1 With Baby Handle
STD-106-12	3 Blades Set 0 ,1 ,2 With Adult Handle
STD-106-13	3 Blades Set 1 ,2 ,3 With Adult Handle
STD-106-14	3 Blades Set 2 ,3 ,4 With Adult Handle
STD-106-15	3 Blades Set 3 ,4 ,5 With Adult Handle
STD-106-16	4 Blades Set 0 ,1 ,2 , 3 With Adult Handle
STD-106-17	4 Blades Set 1 ,2 ,3 , 4 With Adult Handle
STD-106-18	4 Blades Set 2 ,3 ,4 , 5 With Adult Handle
STD-106-19	5 Blades Set 0 ,1 ,2 ,3 , 4 With Adult Handle
STD-106-20	5 Blades Set 1 ,2 ,3 ,4 , 5 With Adult Handle

Single Blades Details			
Code	Blade	B/No	Size
STD:106-00	Premature	00	68 MM
STD:106-0	Neonate	0	80 MM
STD:106-1	Infant	1	95 MM
STD:106-2	Child	2	105 MM
STD:106-3	Medium Adult	3	135 MM
STD:106-4	Large Adult	4	155 MM
STD:106-5	X-Large	5	175 MM

Specification

- SMEGA Hi-light Fiber Optic bundle with 3.5mm dia for improved light transmission and longer life.
- SMEGA LED Technology: For up to 40% more light than conventional halgen bulb
- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- These blades need no maintenance ever
- Both ends of the fiber bundle are protected for life
- 2.5V battery power for even more light
- Compatible with all handles to the ISO 7376 (green Standard)
- Autoclavable.
- Fits all standard F.O. handles (green standard)

Require: 2 x C size batteries Supply set with hard case packing.





LED BLUB

Code :STD 107

Size Measurement

Code	Blade	B/No	Size
STD:107-00	Premature	00	68 MM
STD:107-0	Neonate	0	80 MM
STD:107-1	Infant	1	95 MM
STD:107-2	Child	2	105 MM
STD:107-3	Medium Adult	3	135 MM
STD:107-4	Large Adult	4	155 MM
STD:107-5	X-Large	5	175 MM

Specification

- Innovative Fiber Optic Technology for advanced performance
- SMEGA Hi-light Fiber Optic bundle with 4.3mm dia. And up to 6500 individual micro. Fiber for improved light transmission and longer life.
- SMEGA LED Technology: For up to 40% more light than conventional halgen bulb
- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing. :
- These blades need no 1naintenace ever
- Both ends of the fiber bundle are protected - for life
- 2.5V battery power for even more light
- Compatible with all handles to the ISO 7376 (green Standard)
- Autoclavable.
- Fits all standard F.O. handles (green standard)

Require: 2 x C size batteries Supply set with hard case packing.



SMEGA Miller + Fiber Optic (F.O.) Laryngoscope Blade



LED BULB

Code :STD 108

Size Measurement

Code	Blade	B/No	Size
STD:108-00	Premature	00	65 MM
STD:108-0	Neonate	0	78 MM
STD:108-1	Infant	1	102 MM
STD:108-2	Child	2	154 MM
STD:108-3	Medium Adult	3	195 MM
STD:108-4	Large Adult	4	205 MM

Specification

- SMEGA Hi-light Fiber Optic bundle with 3.5mm dia for improved light transmission and longer life.
- SMEGA LED Technology: For up to 40% more light than conventional halogen bulb
- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- These blades need no maintenance ever
- Both ends of the fiber bundle are protected for life
- 2.5V battery power for even more light
- Compatible with all handles to the ISO 7376 (green Standard)
- Autoclavable.
- Fit all standard F.O. handles (green standard)

Require: 2 x C size batteries Supply set with hard case packing.



SMEGA Macintosh Conventional (C.O.) Laryngoscope Blade



**HALF WELD PIPE
FULL INELD PIPE**



Code :STD 109

Size Measurement

Code	Blade	B/No	Size
STD:109-00	Premature	00	68 MM
STD:109-0	Neonate	0	80 MM
STD:109-1	Infant	1	95 MM
STD:109-2	Child	2	105 MM
STD:109-3	Medium Adult	3	135 MM
STD:109-4	Large Adult	4	155 MM
STD:109-5	X-Large	5	175 MM

Specification

Halogen and LED illumination. Manufactured from AISI 304 stainless steel matt type
Halogen and LED illumination bulb 2.5V apply on the blade easily replace with
few seconds.


C size 2 X Battery handle

: . Compatible with all handles and blades to the ISO 7376 (Standard)

: Autoclavable.

: Fits all standard Conventional (C.O) Handles

Require: 2 x C size batteries Supply set
with hard case packing.





Code :STD 110



LED BULB



Size Measurement

Code	Blade	B/No	Size
STD:110-00	Premature	00	65 MM
STD:110-0	Neonate	0	78 MM
STD:110-1	Infant	1	102 MM
STD:110-2	Child	2	154 MM
STD:110-3	Medium Adult	3	195 MM
STD:110-4	Large Adult	4	205 MM

Specification

- Halgen and LED illumination. Manufactured from AISI 304 stainless steel matt type
- Halgen and LED illumination bulb 2.5V apply on the blade easily replace with few seconds.
- C size 2 X Battery handle
- Compatible with all handles and blades to the ISO 7376 (Standard)
- Autoclavable.
- Fits all standard Conventional (C.O) Handles

Require: 2 x C size batteries Supply set with hard case packing.





LED BULB

Code :STD 111

Size Measurement

Code	Blade	B/No	Size
STD:111-0	Neonate	0	80 MM
STD:111-1	Infant	1	95 MM
STD:111-2	Child	2	105 MM
STD:111-3	Medium Adult	3	135 MM
STD:111-4	Large Adult	4	155 MM

Available Colors



Specification

These disposable stainless steel laryngoscope blades are manufactured for maximum strength and durability. The blades consist of a stainless steel blade and flange with an inter molded durable plastic locking mechanism. The light source utilizes a crystal-clear plastic tubular fiber optic bundle style designed to reflect bright white light from the handle to distal tip or the blade. Compatible with all Green light Fiber Optic systems. Now you can have the durability and strength of a regular stainless steel reusable laryngoscope blade
Compatible with all handles to the ISO 7376 (green Standard)

Require: 2 x C size batteries
Supply set with hard case packing.





LED BULB

Code :STD 112



Size Measurement

Code	Blade	B/No	Size
STD:112-00	Premature	00	68 MM
STD:112-0	Neonate	0	80 MM
STD:112-1	Infant	1	95 MM
STD:112-2	Child	2	105 MM
STD:112-3	Medium Adult	3	135 MM
STD:112-4	Large Adult	4	155 MM

Available Colors



Specification

These disposable stainless steel laryngoscope blades are manufactured for maximum strength and durability. The blades consist of a stainless steel blade and flange with an inter molded durable plastic locking mechanism. The light source utilizes a crystal-clear plastic tubular fiber optic bundle style designed to reflect bright white light from the handle to distal tip or the blade. Compatible with all Green light Fiber Optic systems. Now you can have the durability and strength of a regular stainless steel reusable laryngoscope blade. Compatible with all handles to the ISO 7376 (green Standard)

Require: 2 x C size batteries
Supply set with hard case packing.





Small Medium Stubby

F.O Handle STD 121



Small Medium Stubby

C.O Handle STD 122



**Zipper Bottle
Laryngoscope Packing
Code STD:151**



**Zipper Box
Laryngoscope Packing
Code STD:152**



**Hand Case
Laryngoscope Packing
Code STD:153**



MACINTOSH CONVENTIONAL SETS

Standard Sets

Content	Code
Set Of 3 Blades Size 0, 1, 2 With Adult Handle	STD-170-1
Set Of 3 Blades Size 1, 2, 3 With Adult Handle	STD-170-2
Set Of 3 Blades Size 2, 3, 4 With Adult Handle	STD-170-3
Set Of 3 Blades Size 3, 4, 5 With Adult Handle	STD-170-4
Set Of 4 Blades Size 0, 1, 2, 3 With Adult Handle	STD-170-5
Set Of 4 Blades Size 1, 2, 3, 4 With Adult Handle	STD-170-6
Set Of 4 Blades Size 2, 3, 4, 5 With Adult Handle	STD-170-7
Set Of 5 Blades Size 0, 1, 2, 3, 4 With Adult Handle	STD-170-8
Set Of 5 Blades Size 1, 2, 3, 4, 5 With Adult Handle	STD-170-9
Set Of 6 Blades Size 0, 1, 2, 3, 4, 5 With Adult Handle	STD-170-10
Set Of 6 Blades Size 00, 0, 1 Miller & 2,3,4 Macintosh With Adult Handle	STD-170-11



MILLER CONVENTIONAL SETS

Standard Sets

Content	Code
Set Of 3 Blades Size 00, 0, 1 With Child Handle	STD-170-21
Set Of 3 Blades Size 0, 1, 2 With Child Handle	STD-170-22
Set Of 3 Blades Size 1, 2, 3 With Child Handle	STD-170-23
Set Of 4 Blades Size 00, 0, 1, 2 With Child Handle	STD-170-24
Set Of 4 Blades Size 0, 1, 2, 3 With Child Handle	STD-170-25
Set Of 4 Blades Size 1, 2, 3, 4 With Adult Handle	STD-170-26
Set Of 5 Blades Size 00, 0, 1, 2, 3 With Child Handle	STD-170-27
Set Of 5 Blades Size 0, 1, 2, 3, 4 With Adult Handle	STD-170-28
Set Of 5 Blades Size 1, 2, 3, 4, 5 With Adult Handle	STD-170-29



MCCOY CONVENTIONAL SETS

Standard Sets

Content	Code
Set of 1 Blades Size 2 and Adult Handle	STD-170-31
Set of 2 Blades Size 2,3 and Adult Handle	STD-170-32
Set of 3 Blades Size 1, 2, 3 and Adult Handle	STD-170-33
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	STD-170-34



Set Available In Hard Case Or Soft Leather Pouch
Handle Option Available
Custom Sets Available



MACINTOSH FIBER OPTIC SETS

Standard Sets

Content	Code
Set Of 3 Blades Size 00,0,1 With Baby Handle	STD-180-0
Set Of 3 Blades Size 0, 1, 2 With Adult Handle	STD-180-1
Set Of 3 Blades Size 1, 2, 3 With Adult Handle	STD-180-2
Set Of 3 Blades Size 2, 3, 4 With Adult Handle	STD-180-3
Set Of 3 Blades Size 3, 4, 5 With Adult Handle	STD-180-4
Set Of 4 Blades Size 0, 1, 2, 3 With Adult Handle	STD-180-5
Set Of 4 Blades Size 1, 2, 3, 4 With Adult Handle	STD-180-6
Set Of 4 Blades Size 2, 3, 4, 5 With Adult Handle	STD-180-7
Set Of 5 Blades Size 0, 1, 2, 3, 4 With Adult Handle	STD-180-8
Set Of 5 Blades Size 1, 2, 3, 4, 5 With Adult Handle	STD-180-9
Set Of 6 Blades Size 0, 1, 2, 3, 4, 5 With Adult Handle	STD-180-10
Set Of 6 Blades Size 00, 0, 1 Miller & 2,3,4 Macintosh	STD-180-11



MILLER FIBER OPTIC SETS

Standard Sets

Content	Code
Set Of 3 Blades Size 00, 0, 1 With Child Handle	STD-180-21
Set Of 3 Blades Size 0, 1, 2 With Child Handle	STD-180-22
Set Of 3 Blades Size 1, 2, 3 With Child Handle	STD-180-23
Set Of 4 Blades Size 00, 0, 1, 2 With Child Handle	STD-180-24
Set Of 4 Blades Size 0, 1, 2, 3 With Child Handle	STD-180-25
Set Of 4 Blades Size 1, 2, 3, 4 With Adult Handle	STD-180-26
Set Of 5 Blades Size 00, 0, 1, 2, 3 With Child Handle	STD-180-27
Set Of 5 Blades Size 0, 1, 2, 3, 4 With Adult Handle	STD-180-28
Set Of 5 Blades Size 1, 2, 3, 4, 5 With Adult Handle	STD-180-29



MEGA LIGHT FIBER OPTIC BLADES SET

Standard Sets

Content	Code
Set of 3 Blades Size 00,0,1 and Baby Handle	STD-180-40
Set of 3 Blades Size 0, 1, 2 and Adult Handle	STD-180-41
Set of 3 Blades Size 1, 2,3 and Adult Handle	STD-180-42
Set of 4 Blades Size 0,1, 2, 3 and Adult Handle	STD-180-43
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	STD-180-44



Set Available In Hard Case Or Soft Leather Pouch
Handle Option Available
Custom Sets Available



FO MCCOY BOX MEGA LIGHT BLADE SETS

Standard Sets

Content	Code
Set of 1 Blades Size 2 and Adult Handle	STD-180-61
Set of 2 Blades Size 2, 3 and Adult Handle	STD-180-62
Set of 3 Blades Size 1, 2, 3 and Adult Handle	STD-180-63
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	STD-180-64



FO MCCOY BLADES SETS

Standard Sets

Content	Code
Set of 1 Blades Size 2 and Adult Handle	STD-180-71
Set of 2 Blades Size 2, 3 and Adult Handle	STD-180-72
Set of 3 Blades Size 1, 2, 3 and Adult Handle	STD-180-73
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	STD-180-74



Set Available In Hard Case Or Soft Leather Pouch
Handle Option Available
Custom Sets Available





Section 2

Otoscope

Ophthalmoscope



SMEGA Mini Otoscope SMEGA Fiber Optic (F.O) Mini Otoscope



NEW

Code :STD 201



NEW

Code :STD 202



Available Colors



Specification

Modern pocket otoscope in a unique, compact design. Handy and powerful. Instrument head matt-black inside. Eliminates reflexes. Distal Fiber-Optic (F.O) Illumination with Enhanced LED Technology with 100% more light compared with, LED lamp. Attachment clip with integrated on/off switch. Secure. Switches off automatically when returned to the pocket

Supply: 10 single use mini speculum (Ø2.5 and 4mm) with storage bottle packing in black pouch.
Require: 2x AA batteries (not included with set)

Specification

Modern pocket otoscope compact design. Handy and powerful. Distal Fiber-Optic (F.O) illumination 4mm bundle with Enhanced LED Technology with 100% more light compared with LED Lamp. Attachment clip with integrated on/off switch.

Supply: 10 single use mini speculum (Ø2.5 and 4mm) with storage bottle packing in black pouch.
Require: 2x AA batteries (not included with set)



SMEGA Mini Otoscope SMEGA Fiber Optic (F.O) Mini Otoscope



Specification

High quality heine type plastic Otoscope with extra bright light 2.5v LED Illumination bulb illumination provided thought fiber optic system. High quality swiveling viewing window with 3x magnification lens, out put for pneumatic membrane test. Easily insert battery facility in handle

Supply: 10 single use mini speculum (2.5 and 4mm) with storage bottle packing in black pouch.
Require: 2x AA batteries (not included with set)

Specification

High quality professional plastic Otoscope with extra bright light 2.5v LED Illumination bulb illumination provided thought fiber optic system. High quality swiveling viewing window with 3x magnification lens, out put for pneumatic membrane test.

Supply: 10 single use mini speculum (Ø2.5 and 4mm) with storage bottle packing in black pouch.
Require: 2x AA batteries (not included with set)



SMEGA Otoscope Fiber Optic For Adult & Child



NEW

Code :STD 205



NEW

Code :STD 206



Available Colors



For Adult

Specification

- Bayonet fitting
- Illumination light 2.5v LED bulb
- Metal handle C
- Plastic Head CPP
- Rheostat dimmer
- 05 x Ø2.5mm
- 05 x Ø4.0mm

Supply: in zipper pouch packing
Require: 2 x C batteries (not included with Set)

For Child

Specification

- Click fitting
- illumination light 2.5v LED bulb
- Metal battery child chrome plated
- Plastic Head CPP
- Reusable handle C
- Dimmable rheostat
- Disposable ear funnel
- 05 x Ø2.5mm
- 05 x Ø4.0mm

Supply: in zipper pouch packing
Require: 2 x AA batteries (not included with Set)



SMEGA Fiber Optic (F.O) Adult Otoscope



Code :STD 207

Code :STD 208

Available Colors



Specification

High quality professional plastic head, chrome plated battery handle Otoscope with extra bright light 2.5v LED bulb 2.5v illumination provided thought fiber optic system. High quality Big Swiveling viewing window with 3x magnification lens output for pneumatic membrane test,

Supply: 10 single use mini speculum (Ø 2.5 and 4mm) with storage bottle packing in black pouch.
Require: 2x C batteries (not included with, set)

Specification

Parker F.O. Otoscope set high quality professional metal head, chrome plated battery handle Otoscope with extra bright light 2.5v LED bulb illumination provided thought fiber optic system. High quality big swiveling viewing window, with 3x magnification lens output for pneumatic mc1nbrane test.

Supply: 10 single use mini speculum (0 2.5 and 4mm) with storage bottle packing in black pouch.
Require: 2x C batteries (not included with, set)



SMEGA Conventional (C.O) Otoscope



Code :STD 209

Code :STD 210

Available Colors



Specification

SMEGA PARKER OTOSCOPE

Realizable and indestructible high quality instruments Invaluable for everyday ear throat and nose use Regulated easily neutral 2.5v LED white light
3 speculum (size Ø 2.5 3.5 and 4.5mm)
3x High quality magnifying lens

Available: Polish and Sand finish
Require: 2x C batteries (not included with set)

Specification

SMEGA GM Otoscope Bayonet connecting system head fitting ideal for pneumatic tympanic membrane test can also be used for rhinoscopic examinations.
2.5v LED white light
3 speculum (size Ø 2.5 3.5 and 4.5mm)
3x High quality magnifying lens

Available: Polish and Sand finish
Require: 2x C batteries (not included with set)



SMEGA Conventional (C.O) Otoscope



NEW

Code :STD 211



Code :STD 212



Available Colors



Specification

- Screw fitting
- illumination light 2.5v bulb LED
- Metal handle C
- Plastic Head CPP
- Rheostat dimmer
- 05x Ø 2.5mm
- 05x Ø 4.0mm
- Supply: in pouch packing

Specification

2.5v standard illumination LED bulb. Chrome plated handle and head 4x magnifying glasses. Reusable speculum size Ø2.5 3.5 and 4mm include. Connecting for pneumatic Otoscope. Battery handle with rheostat for adjustment of illumination.

Supply: German bulb Technology for super view with 3 reusable ear specula.

Packing: Soft pouch and plastic hard case
Require: 2x C batteries (not included with set)
Available: Polish and Sand finish



SMEGA Conventional (C.O) Otoscope



NEW

Code :STD 213



Code :STD 214



Available Colors



Specification

Modern pocket otoscope in unique, compact design. compared with conventional bulbs. Viewing, window with 3X Switches off automatically when returned to the pocket.

2-piece handle and head. Easy to maintain, flexible

Supply: 10 single use mini speculum Ø2.5 and 4mm with storage bottle packing in black pouch.

Require: 2x AA batteries (not included with set)

Specification

Complete plastic Otosoope with 3 specula size Ø2.5 3.5 and 4.5mm

2.5v LED light bulb

Rheostat adjustable light

Big viewing window with 3x magnification lens

Packing: Soft pouch and plastic hard case

Require: 2x C batteries (Not included with set)



SMEGA Mini Ophthalmoscope



Specification STD 219

Top quality ophthalmoscope with fixation star aperture in pocket format. Quality optics in mini-format. Ideal for rounds and house calls. Attachment clip with integrated on/off switch. Secure. Switches off automatically when returned to the pocket. High-quality handle: White upper section/ refined plastic. Shockproof, sturdy, non-slip. Enhanced LED Technology. Concentrated white light. for perfect illumination and a brilliant image. 4 different apertures. A full-function instrument incl. fixation star. Range of lenses: -20 D to +20 D. 2-piece handle and head. Easy to maintain, flexible. Replaceable batteries. Size AA
Supply: Bristle and Hard plastic case packing,

Specification STD 220

(small, large, blue, half moon focus) 2.5v LED and LED lamp light for true tissue color and consistent long lasting illumination 16 color coded lenses -25 to +25 for excellent resolution (Red for - black for+) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification of adopter selling.

Supply: Bristle and hard plastic case packing,
Require: 2x AA batteries (not included with set)



SMEGA Mini & Standard Ophthalmoscope



Code :STD 219



Code :STD 220

Specification

GM Ophthalmoscope set pin contacting fitting
Ophthalmoscope head: pin contact fitting. Lens selection wheel permits choice or corrective lenses from 0 to 25 and 0 to + 40 diopters single large circle aperture. Rheostat for adjustment of illumination 2 battery handle choice

Require: 2 x C size batteries (not included w/U1 set)
Supply: Hard case.

Specification

Parker Ophthalmoscope 2.5v LED bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses. -25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear Identification fo adopter setting.

Supply: Soft pouch or plastic hard case
Require : 2 x C batteries (not included with set)

