



EN Surgical Blades & Sutures

DE Chirurgische Klingen & Nähte

ES Cuchilla & Suturas quirúrgica

Page 801 to 803




- Item Codes for disposable surgical blades are as follows:
HY-DSB10 through HY-DSB36d plus SS or CS, respectively (eg.: HY-DSB13-SS)



Item Codes for sterile stitch cutters are:
HY-SSC065-CS, HY-SSC110-CS, HY-SMTB

For the Item Codes of sterile disposable scalpels
choose your blade type above: HY-SDS13-SS



Item Codes are based on material and if there is a needle or not:

HY-PGA-N with needle, HY-PGA-O without needle (please check page 803 for needle types)

ABSORBABLE SUTURES



Polyglycolic Acid (PGA)

Synthetic, Absorbable, Braided Suture
With Needle or Without Needle



Poliglecaprone 25 (MO)

Synthetic, Absorbable, Monofilament Suture
With Needle or Without Needle



Rapid Polyglycolic Acid (PGAR)

Synthetic, Absorbable, Braided Suture
With Needle or Without Needle



Chromic Catgut (CC)

Absorbable Suture
With Needle or Without Needle



Polyglactin 910 (PGLA)

Synthetic, Absorbable, Braided Suture
With Needle or Without Needle



Plain Catgut (CP)

Absorbable Suture
With Needle or Without Needle



Polydioxanone (PDS)

Synthetic, Absorbable, Monofilament Suture
With Needle or Without Needle

NON-ABSORBABLE SUTURES



Silk (SK)

Non-absorbable, Braided, Natural Suture
With Needle or Without Needle




Polyester (PB)

Synthetic, Non-absorbable, Braided Suture
With Needle or Without Needle



Nylon (NL)

Synthetic, Non-absorbable, Monofilament
Suture With Needle or Without Needle




Stainless Steel Wire (SS)

Synthetic, Non-absorbable, Monofilament
Suture With Needle or Without Needle




Polypropylene (PM)

Synthetic, Non-absorbable, Monofilament
Suture With Needle or Without Needle




Item Codes for needles with Basic Shape and Major Type codes are added to the suture item code above:
HY-PGA-1/2C-T1 or HY-PGA-SKI-T1 (indicate needle length, diameter, etc. if needed)

Needle Anatomy



Basic Shapes



Major Types

Needle Shape (针形)	Point Type	Symbol (符号)
	Round Bodied	● T1
	Curved Cutting	▲ T2
	Reverse Cutting	▼ T3
	Tapercut	⊙ T4
	Micro-point Reverse Cutting	▽ T5
	Micro-point Spatula Curved	◀ T6
	Blunt Taper Point	⊙ T7