

ECRS

OPHTHALMOLOGY
CATARACT
REMOVAL SET
GA8363521-55129



2
0
2
1

Hygeia Medical



CODE	INSTRUMENT	REMARKS	SIZE
TCF-04-13B	Castroviejo Needle Holder	with lock	(12.5cm)
TCF-04-13	Barraquer Needle Holder		(15cm)
DF-371A-4331A	Castroviejo Scissors	left	(10cm)
DF-371A-4331B	Castroviejo Scissors	right	(10cm)
DF-371B-4315M	Castroviejo Corneal Scissors		(10cm)
DF-371B-4316P	Aebli Capsule Scissors		(10cm)
DF-4-132	Scalpel Handle No.3	solid	(12 cm)
DF-360-4214	Graeffe Forceps	1x2 teeth	(7cm)
DF-362-4235	Gill-Hess Iris Forceps		(9cm)
DF-362-4234	Gill Iris Forceps		(9cm)
DF-362-4242	Bishop-Harmann Fixation Forceps	1x2 teeth	(8.5cm)
DF-365-4266	Castroviejo Suture Forceps		(10cm)
DF-366-4278	Bonaccolto Suture Tying Platform	smooth, angled	(11cm)
DF-362-4239	Hess Forceps	double angled	(7cm)
2018-EOL-M	Kelman-McPherson Tying Platform	angled	
2019-EOL-M	Masked Capsorheads	very delicate sharp tips	
2020-EOL-M	Dardenne Suture Tying Platform	6mm angeled	
DF-341-3892	Iris Knife	2x5mm cutting edges	(12cm)
DF-348-4020	Stevens Hook		(12cm)
DF-351-4074	Lewis Lens Scoop	smooth	(12cm)
DF-350-4059	Barraquer Iris Spatula	smooth	(12cm)
DF-350-4053	Castroviejo Spatula		(12cm)
DF-350-4049	Jaffe Lens Spatula		(12cm)
DF-350-4050	Kwitko Spatula		(12cm)
DF-350-4048	Hirschmann Spatula		(12cm)
DF-350-4065	Castroviejo C Spatula		(10cm)
DF-354-4127	Wilder		(10cm)
DF-333-3778	Barraquer Colibri	(adult eye lid)	(4cm)
DF-333-3777	Barraquer Colibri	(child size lid)	(3cm)
FSI-AIC-27	Air Injection Cannula	angled 5mm at blunt tip	27 gauges
FSI-RCC-21	Randolph Cyklodialyse Cannula	angled 13mm f. flat tip	21 gauges
FSI-LC-23	Lacrimal Cannula	malleable tip	23 gauges
FSI-SIC-23	Simcoe Irrigation Cannula		23 gauges
DF-337-3818	Desmarres Lid Retractor		(18cm)



Hygeia Medical Ltd. (Eire)
KYLETELOGUE HOUSE - 9 CHERRYWOOD – TALLANSTOWN
DUNDALK 21222 - CO. LOUTH – IRELAND

Hygeia Instruments (Malaysia) Sdn Bhd
3A-4-139 KOMPLEX BUKIT JAMBUL – JALAN RUMBIA
11900 BAYAN LEPAS – PU:AU PINANG – MALAYSIA

www.hygeiainstruments.icu



TCF-04-13B
Castroviejo
Needle Holder
with lock
(12.5cm)



TCF-04-13
Barraquer
Needle Holder
(15cm)



DF-371A-4331A
Castroviejo Scissors
left (10cm)



DF-371A-4331B
Castroviejo Scissors
right (10cm)



DF-371B-4315M
Castroviejo
Corneal Scissors
(10cm)



DF-371B-4316P
Aebli Capsule Scissors
(10cm)



DF-4-132 Scalpel Handle No. 3 – solid (12cm)



DF-360-4214
Graeffe Forceps
1x2 (7cm)



DF-362-4242
Bishop-Harmann
Fixation Forceps
1x2 (8.5cm)



DF-365-4266
Castroviejo Suture Forceps
(10cm)



DF-362-4235
Gill-Hess Iris Forceps
(9cm)



DF-362-4234
Gill Iris Forceps
(9cm)

2018-EOL-M
Kelman-McPherson
angled tying platform

2019-EOL-M
Masket Capsulorheads
very delicate sharp tips

2020-EOL-M
Dardenne Suture Tying
6mm tying platform
(10cm)

DF-366-4278
Bonaccolto
Suture Tying Platform
smooth, angled
(11cm)

DF-362-4239
Hess Forceps
(7cm)





DF-341-3892
Iris knife (12cm)
with 2 cutting edge
of 5mm



DF-348-4020
Stevens Hook
(12cm)




DF-351-4074
Lewis Lens Scoop
smooth (12cm)



DF-350-4059
Barraquer Iris Spatula
smooth (12cm)



DF-350-4053
Castroviejo Spatula
(12cm)



DF-350-4049
Jaffe Lens Spatula
(12cm)



DF-350-4050
Kwitko Spatula
(12cm)



DF-350-4048
Hirschmann Spatula
(12cm)



DF-350-4065
Castroviejo C Spatula
(10cm)



DF-354-4127
Wilder
(10cm)



DF-333-3778
Barraquer
Adult size (4cm)

DF-333-3777
Barraquer (Colibri)
Child size (3cm)



DF-337-3818
Desmarres
Lid Retractor
(18cm)



FSI-AIC-27
Air Injection Canula
angled 5mm at blunt tip
27 gauges



FIS-RCC-21
Randolph Cyklodialyse Cannula
angled 13mm from flat tip
21 gauges



FSI-LC-23
Lacrimal Cannula
malleable tip
23 gauges



FSI-SIC-23
Simcoe Irrigation Cannula
23 gauges